

# KU-2™ Viscometer

for Paints, Coatings, and Inks

**ASTM D562 Compatible**  
(industry specification)

**Easy to use**  
no weights, simplifies an established test procedure

**LED Display Info:**

- Krebs Units
- Gram Units (Weight)
- Centipoise\*

Select Krebs or Grams or Centipoise

**Lock-In Test Results**  
with Hold Switch

**Accuracy: ±1.0% of range**

**Repeatability: ±0.5%**

**Standard Krebs Spindle**

**Measurement range:**  
40 to 141 KU, 32 to 1099 gm,  
and 27 to 5274 cP\*

**Printer Connection**  
for automatic test documentation

**Adapter provided**  
for 1/2 Pint, Pint,  
& Quart Containers



## What's Included?

- Instrument
- Krebs-type Spindle (p48)
- Communications Port for Printer
- Adapter for Quart, Pint and Half-Pint Cans†

## Optional Accessories

- Air Purge
- Paste Spindle (p48)

## Applications

- Paints
- Coatings
- Adhesives
- Inks
- Pastes

\*Centipoise values based on the conversion from Krebs Units as defined in the ASTM D562.

†Adapters which accommodate metric size cans are available for an additional fee — contact Brookfield to discuss your requirements